



Panorama™
Next-generation NIPT

A pregnant woman with her hair in a bun, wearing a bright yellow two-piece outfit, stands in a field. A complex digital network of white lines and blue and green 'x' markers is overlaid on the image, resembling a neural network or data flow.

See more with Panorama AI

The most widely used NIPT is now
powered by Artificial Intelligence

The only SNP-based NIPT delivers more insights and greater accuracy

Panorama's single nucleotide polymorphism (SNP)-based non-invasive prenatal test (NIPT) is:


- the most rigorously validated NIPT
- the only NIPT that distinguishes mother's DNA from baby's DNA
- the approach that creates unique, clinically validated capabilities

The #1 NIPT is now powered by Artificial Intelligence


Panorama AI leverages artificial intelligence (AI) to learn from the more than 2 million tests already processed by Natera.

Panorama AI combines AI with Natera's proprietary SNP-based methodology to give accurate results on difficult-to-call cases.

- Maintaining industry-leading accuracy while significantly lowering "no-call" rates
- Increasing accuracy for 22q11.2 deletion detection¹



SNPs
Panorama evaluates SNPs — the 1% of our DNA that makes us different from one another.



#1 for a reason¹³

TESTS
2M+
performed

PATIENTS
1.3M+
studied

SNP TARGETS
13,392

PUBLICATIONS
23+




Panorama AI was validated in SMART, the largest prospective NIPT study^{1,5,6}



“No-calls” provide actionable information

SMART showed that, for Panorama, patients with “no-call” results after redraw have significantly higher rates of adverse outcomes—compared to baseline or patients with a call after redraw.⁶

Adverse pregnancy outcome rates for Panorama (prior version)



CLINICAL UTILITY OF “NO-CALLS”

Panorama AI further enriches this increased-risk group, making the “no-call” after redraw group highly actionable.^{6,13}

Consider the following for this group:

- Increased surveillance for pre-term labor and preeclampsia



“Patients with failed cfDNA should have follow-up surveillance.”

MARY NORTON, MD

Based upon data presented by Dr. Mary Norton at SMFM 2021.⁶
*Includes intrauterine fetal demise, stillbirth, miscarriage, and termination.

Highest accuracy. Lower “no-calls.”


Historically, there has been a trade-off between NIPT accuracy and “no-call” rates. With Panorama AI, that paradigm has been broken, delivering the best of both worlds.^{1-4,14,15,18-20,23-25}



Unparalleled screening for 22q11.2

Panorama AI has increased accuracy for 22q11.2 deletion screening.⁵

Panorama AI can identify most 22q11.2 deletions, including small deletions of 0.5 Mb.⁵ Most NIPTs can only screen for large microdeletions—above 2.5 Mb.^{22,27}




>99.9% Sensitivity⁵ Full A–D deletion (2.5 Mb+)

83% Sensitivity⁵ Full & nested deletions (0.5 Mb+)

53% PPV⁵

22q11.2 had an incidence rate that was comparable to some of the common trisomies^{1,5}




ACTIONABLE PERFORMANCE

Panorama AI has increased sensitivity and >2X the positive predictive value (PPV) for 22q11.2 deletions, enhancing the impact of this screen, given the number of prenatal/neonatal interventions for this disorder.⁵


*Most other studies of this kind rely on observational assessment, not genetic diagnostic testing.

Panorama maintains industry-leading performance and clinical value in singleton pregnancies

Combined FNR in clinical validation studies (T21, T18, T13)



Fetal sex error rate in clinical validation studies



Panorama provides unique information to appropriately triage high-risk twin pregnancies

Zygoty is key in appropriately managing twin pregnancies and monitoring for such conditions as twin-twin transfusion syndrome (TTTS):

- One-fifth of monozygotic twins are misidentified by ultrasound.²⁸
- One-sixth of misidentified cases that are referred to fetoscopic centers reach stage IV TTTS, and one-fifth suffer fetal/neonatal demise.²⁹

According to ACOG and SMFM, divergent fetal fraction in dizygotic twins can make it difficult to detect aneuploidy in a twin with low fetal fraction if the other twin is euploid and has a high fetal fraction.³⁰

Only Panorama NIPT can detect zygoty

ZERO errors in zygoty detection^{7*}

Only Panorama reports individual fetal fraction**

7% of dizygotic twins have highly divergent fetal fraction⁸

Every pregnancy deserves Panorama

“(NIPT is) the most sensitive and specific screening test for the common aneuploidies...should be discussed and offered to all patients regardless of maternal age or baseline risk.”³⁰

ACOG PRACTICE BULLETIN 226

ACOG guidelines highlights, including unique capabilities of SNP-based NIPT

ACOG PRACTICE BULLETIN 226 HIGHLIGHTS ³⁰	PANORAMA SNP-BASED NIPT ^{1-4,7-12}	COUNTING-BASED NIPT METHODS ^{14,15,18-20,23-25}
“Cell-free DNA [NIPT] can be performed in twin pregnancies”	✓	✓
“...one laboratory method which uses SNP analysis reports zygoty as well as individual fetal fractions.”	✓	✗
“[In twin gestations,] it is possible that an aneuploid fetus would contribute less fetal DNA , therefore masking the aneuploid result.”	✓	✗
Highlights known sources of false positives for most NIPTs: • Vanishing twin • Maternal mosaicism	Addressed	Not addressed
“Of the [NIPT] methods, the [SNP-based] method can identify triploidy.”	✓	✗
“...the lower prevalence of fetal trisomies in younger patients results in...a lower PPV...”	Robust clinical PPVs published	Analytical PPVs available


Panorama's areas of unique clinical value

Twins differentiation

Every result matters

Accuracy and reliability

Trisomy 21 PPVs by age groups, shown in a large, 1M patient outcomes study^{12,13,31,32}



PPVs

Strong clinical evidence with robust PPVs, regardless of age.

*In validation studies
**For dizygotic twins



NateraCore

— simple, tailored resources to support you and your patients every step of the way



<p>Education Patient-friendly materials and information sessions, covering basic genetics to specific tests</p>	<p>Access Programs and price transparency – rooted in our commitment to provide affordable testing for all who can benefit</p>	<p>Ordering Flexible options based around your needs, including intuitive remote ordering</p>	<p>Results Clear, actionable reports, served with time-saving tools and a side of expert guidance</p>	<p>Next steps Value-add services that go beyond the test to address what's next</p>
--	---	--	--	--

A selection of our NateraCore offerings

<p>Pre- and post-test genetic information sessions – access to board-certified genetic counselors, available to all providers and patients</p> 	<p>Price Transparency Program (PTP) – personalized cost estimates and a self-pay cash alternative</p> 
<p>Virtual testing – fully remote testing option, combining online ordering with self-service patient education and mobile phlebotomy</p> 	<p>Flexible phlebotomy options – via local, Natera-approved lab or at-home mobile services; available in all states, at no cost to patients</p> 

Panorama screens for:

Singleton pregnancies

- Trisomies 21, 18, 13
- Monosomy X
- Triploidy
- Sex chromosome trisomies*
- 22q11.2 deletion syndrome (optional)
- Additional microdeletion syndromes (optional)
- Fetal sex (optional)

Twin pregnancies

- Zygosity
- Trisomies 21, 18, 13
- Fetal sex for each twin (optional)

If screening reveals monozygotic twins, Panorama can additionally screen for:

- Monosomy X
- Sex chromosome trisomies*
- 22q11.2 deletion syndrome (optional)

Egg-donor or surrogate pregnancies (singleton pregnancies only)

- Trisomies 21, 18, 13
- Fetal sex (optional)

References

- Dar et al. Multicenter prospective study of SNP-based cfDNA screening for aneuploidy with genetic confirmation in 18,497 pregnancies. *Society of Maternal-Fetal Medicine, SMFM*. Virtual Meeting. Oral Presentation. Jan 25-30, 2021
- Pergament et al. *Obstet Gynecol*. 2014 Aug; 124(2 Pt 1):210-8
- Nicolaides et al. *Prenat Diagn*. 2013 June; 33(6):575-9
- Riyan et al. *Fetal Diagn Ther*. 2016;40(3):219-223
- Dar et al. Multicenter prospective study of SNP-based cfDNA for 22q11.2 deletion in 18,290 pregnancies with genetic confirmation. *Society of Maternal-Fetal Medicine, SMFM*. Virtual Meeting. Oral Presentation. Jan 25-30, 2021
- Norton et al. Perinatal and genetic outcomes associated with no call cfDNA results in 18,497 pregnancies. *Society of Maternal-Fetal Medicine, SMFM*. Virtual Meeting. Oral Presentation. Jan 25-30, 2021
- Norwitz et al. *J Clin Med*. 2019 Jun; 8:937
- Hedriana H et al. *Prenat Diagn*. 2020 Jan;40(2):179-184
- Nicolaides et al. *Fetal Diagn Ther*. 2014;35(3):212-7
- McKanna T et al. *Ultrasound Obstet Gynecol*. 2019;53(1):73-79
- Martin KA et al. *Am J Obstet Gynecol*. MFM 2020;2:100152
- DiNunno W. et al. *J Clin Med*. 2019 Aug; 8,1311
- Natera internal data on file
- Stokowski et al. *Prenat Diagn*. 2015 Dec; 35(12):1243-6
- Jones et al. *Ultrasound Obstet Gynecol*. 2018 Feb;51(2):275-276
- Hooks et al. *Prenat Diagn*. 2014;34(5):496-499
- Schmid et al. *Fetal Diagn Ther*. 2017, DOI: 10.1159/000484317
- Porreco et al. *Am J Obstet Gynecol*. 2014;211:365.e1-12
- Mazloom et al. *Prenat Diagn*. 2013 Jun;33(6):591-7
- Tynan et al. *Society for Maternal-Fetal Medicine, SMFM*. Las Vegas, Nevada. Jan 23-28, 2017
- Sehnert et al. *Clin Chem*. 2011 Jun;57(7):1042-1049
- Bianchi et al. *Obstet Gynecol*. 2012 May; 119(5):890-901
- Bianchi et al. *N Engl J Med*. 2014;370:799-808
- Verinata white paper. Analytical validation of the Verifi prenatal test. 2012
- Commercial protocol not validated; Illumina marketing materials cite "Srinivasan et al. *Am J Hum Genet*. 2013 Feb 7; 92(2): 167-176" which does not match number of reads used in commercial testing
- Blumenfeld et al. *J Ultrasound Med*. 2014 Dec;33(12):2187-92
- Baud et al. *Ultrasound Obstet Gynecol*. 2014; 44: 205-209
- ACOG Practice Bulletin 226. *Obstet Gynecol*. 2020 Oct;136(4):859-867.
- Myriad Women's Health website accessed Oct 23 2019
- Progenity Innatal Clinician Guide.

13011 McCallen Pass, Building A Suite 100 | Austin, TX 78753 | natera.com

The tests described has been developed and its performance characteristics determined by the CLIA-certified laboratory performing the test. The test has not been cleared or approved by the US Food and Drug Administration (FDA). Although FDA is exercising enforcement discretion of premarket review and other regulations for laboratory-developed tests in the US, certification of the laboratory is required under CLIA to ensure the quality and validity of the tests. CAP accredited, ISO 13485 certified, and CLIA certified.
© 2021 Natera, Inc. All Rights Reserved. PAN_BR_MD_Pano_AI_20210312_NAT-8020378



*Reported when suspected