



FARAPULSE™

Pulsed Field Ablation System






Safety Profile to Adjacent Structures


Thermal Energy Can't Differentiate; PFA is Tissue Selective

Thermal lesions can spread into surrounding tissue indiscriminately

RADIOFREQUENCY ABLATION



CRYO ABLATION




Potential Complications:

- PV Stenosis
- AE Fistula
- Phrenic Nerve Palsy


¹Reddy VY et al. J Am Coll Cardiol. 2019;74(3): 315-26 [LINK](#)

PFA can ablate myocardium while reducing collateral damage¹

PULSED FIELD ABLATION



Cardiomyocytes have low thresholds to PFA




Other tissue/cell types are more resistant to PFA & remained uninjured despite exposure to the field



Pulse Field Ablation Mechanism of Action

Non-Thermal IRE

Irreversible Electroporation: cell membrane pore formation & cell death occurs with sufficient electric field



T. Kotnik et al, *IEEE Electrical Insulation Magazine*, Vol. 28, No. 5 p. 14-23, 2012


Caution: Investigational Devices. Limited by Federal (or US) law to investigational use only. Not available for sale.
Boston Scientific Confidential – Access Limited to Authorized Personnel. Do not copy, display or distribute externally.



FARAPULSE Dedicated History

Pioneering PFA Technology

FARAPULSE has maintained a singular focus on cardiac PFA since 2012



Design Goals:



BUILD

the ideal PFA system from the ground-up



OPTIMIZE

through pre-clinical & clinical safety, efficacy & lesion durability



VALIDATE

system in randomized trials & commercial use



The FARAPULSE™ PFA System

Boston
Scientific

Design Goals

FARADRIVE™*


Steerable Sheath



Steerable sheath designed to navigate the FARAWAVE catheter to targeted cardiac anatomy

FARAWAVE™*

Pulsed Field Ablation Catheter



Over-the-wire PFA catheter with distal end variable morphology **designed to treat a range of PV anatomies**

FARASTAR™*

Pulsed Field Ablation Generator











Bipolar & biphasic waveform generator designed to deliver proprietary PFA therapy through the FARAWAVE catheter



Clinical Data – Lesion Durability at Remap

Acute Isolation Does Not Equal Durable Lesions

110 patients returned for prospective remaps at 93 ± 30 days

Cohort	N	Acute PVI (% PVs)	~3 Month Remap		
			Durable PVI (% PVs)	Durable PVI (% pts)	
 Monophasic General anesthesia & paralytics	11	100%		45%	18%
 Biphasic (Early/Other) Waveform eliminated need for GA	55	100%		84%	58%
 Optimized Biphasic 4 paired applications (8 total) per PV at 1.8 – 2.0kV <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">  <p>2 pairs in 'Basket'</p> </div> <div style="text-align: center;">  <p>2 pairs in 'Flower'</p> </div> </div>	44	100%		96%	84%

FARAPULSE Optimized Biphasic 96% PVs Durable Isolation at ~3 months¹

¹Reddy et al. JACC: Clinical Electrophysiology 7.5 (2021): 614-627 [LINK](#)


Caution: Investigational Devices. Limited by Federal (or US) law to investigational use only. Not available for sale.
 Boston Scientific Confidential – Access Limited to Authorized Personnel. Do not copy, display or distribute externally.



Real-World EU Single Center Experience

Experienced Short Learning Curve & Efficient Procedure

Learning curve of PFA Ablation



Procedural characteristics, n	Phase 1	Phase 2	P value	Overall
	25 Pts, 98 PV	166 Pts, 650 PVs		191 Pts, 748 PVs
Duration of overall procedure, mins				
Mean±SD, n	46±14	38±14	0.004*	39±14
Median (min-max)	43 (28-92)	35 (19-90)		36 (19-92)
Overall fluoroscopy time, mins				
Mean±SD, n	11±5	8±4	0.001*	9±4
Median (min-max)	9.9 (5.1-22.9)	7.5 (3.9-27.9)		7.7 (3.9-27.9)
Number of catheters used				
1, n (%)	24 (96%)	166 (100%)		190 (99.5%)
2, n (%)	1 (4%)	0 (0%)		1 (0.5%)
Catheter size (31/35 mm), n (%)	14/11 (56%/44%)	82/84 (49%/51%)		96/95 (50%/50%)
Number of veins attempted, n	98	650	...	748
Number of applications, per vein	8±1	8±0		8±1
Single-shot isolation,† n (%)	95 (97%)	649 (100%)	0.007‡	744 (99.5%)
Patients with all PVs single-shot isolation, n (%)	22 (88%)	165 (99%)	0.006‡	187 (98%)



Committed to Advancing PFA

Near term clinical studies and form factors

Paroxysmal Indication

FARAVIEW
Integrated Mapping


Persistent Indication &
Linear PFA Catheter

Front-Line Persistent

ADVENT




FARAWAVE



ADVANTAGE


+



FARAPOINT

**AVANT
GUARD**





Harnessing the power of **RHYTHMIA** + **FARAWAVE**

RHYTHMIA HDx Integration Design Goals**

Catheter visualization with reduced fluoroscopy

Workflow support

Advanced versatility beyond PVI



An Integrated Experienced FARAWAVE™ with FARAVIEW™ Mapping Module**

Boston
Scientific

Fluoro-Based / Impedance Tracked




Pulmonary
Vein Isolation

Posterior Wall
Isolation



FARAVIEW™ Mapping Module






Maintain Workflow Simplicity

FARAWAVE™ with FARAVIEW™ Mapping Module**


Boston
Scientific



**Construct Geometry with
FARAWAVE**



**Position FARAWAVE within
the Target Vein**



**Visualize Rotation of the
Electric Field**


*Tag shape for illustration /
workflow guidance only*




Visualize the Electric Field** Field Tags Overview

Field Tags are a new tagging methodology designed specifically for PFA


Approach:



Estimate the
Field Volume






Visualize
where the field
intersects with
the anatomy



Automatically
tag highlighted
area when PFA
is detected

Saturation of Field Tag color increases
with overlapping ablations







RHYTHMIA HDx™ Integration FARAWAVE™ Workflow Guidance**

Boston
Scientific

LA mold on Bench



Workflow guidance through tagging (and tag overlap) to assess uniformity of the ablation set.



Appendix – Registry Data
MANIFEST-PF & EU-PORIA
VersaCross Connect



MANIFEST-PF (Safety and Effectiveness of Pulsed Field Ablation to Treat Atrial Fibrillation: One-Year Outcomes from the MANIFEST-PF Registry)



DOI: 10.1161/CIRCULATIONAHA.123.064959

Turagam, Mohit K., et al. "Safety and Effectiveness of Pulsed Field Ablation to Treat Atrial Fibrillation: One-Year Outcomes from the MANIFEST-PF Registry." *Circulation* (2023).

Objective

Retrospective analysis of safety, procedural endpoints, and long-term outcomes of all consecutive FARAPULSE™ cases, inclusive of learning curve.

Registry Design



24 European Centers



77 Operators



1568 Patients
65% PAF



All FARAPULSE™ PFA
Cases from 3/21-5/22

Results



61 min Median
Procedure Time



1.9% Major Complication Rate
4.0% Minor Complication Rate

81.6%

PAF One-Year
Freedom from AF/AFL/AT

Conclusions

- In an all-comer AF patient population of the first patients undergoing ablation with FARAPULSE
 - The acute PV isolation rate was 99.2% with ~1 hour procedure times
 - The one-year freedom from AF/AFL/AT was 81.6% for paroxysmal AF patients
 - The adverse event rate was low with no reported esophageal damage or PV stenosis

Full Cohort and AF Sub-Type

78.1% Full Cohort

81.6% PAF

71.3% PerAF

73.5% LS PerAF

By Center FARAPULSE™

Procedural Volume

(PFA procedure number)

78.7% >100, n = 802

79.1% 50-100, n = 350

76.0% 0-50, n = 416



EU-PORIA (European Real-World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation)

DOI: 10.1093/europace/euad185

B. Schmidt, et al. "European Real World Outcomes with Pulsed Field Ablation in Patients with Symptomatic Atrial Fibrillation - Lessons from the multicenter EU-PORIA Registry." *Europace* (2023).

Objective

To describe real-world adoption, workflow, acute and long-term outcomes after pulsed field ablation (PFA) in an all-comer atrial fibrillation (AF) patient population in high-volume European centers, inclusive of learning curve.

Registry Design



7 High Volume Centers
(400-1400 ablations/yr)



42 Operators



1233 Patients
60% PAF



All FARAPULSE™ PFA
Cases from 3/21-5/22

Results



58 min Median
Procedure Time



1.7% Major Complication Rate
1.9% Minor Complication Rate

80%

PAF One-Year
Freedom from AF/AT

Table 1 : Freedom from AF/AT* and Vein Reconnection rate

One-Year Freedom from AF/AT	74% [95% CI, 71%-76%], full cohort 80% [95% CI, 77%-83%], paroxysmal AF 66% [95% CI, 61%-71%], persistent AF 67% [95% CI, 48%-82%], LS persistent AF
One-Year Freedom from AF/AT by Operator AF Ablation Experience	65%, [95% CI, 31%-88%] <2 years 72%, [95% CI, 66%-77%] 2-5 years 76%, [95% CI, 73%-79%] >5 years
Freedom from AF/AT Recurrence by Previous Primary Ablation Modality	75%, [95% CI, 70%-80%] Radiofrequency (RF) 72%, [95% CI, 65%-78%] Cryoablation (cryo) 75%, [95% CI, 72%-79%] Both RF and cryo
Vein Reconnection Rate during Repeat Ablation	149/1233 repeat ablations at 226 [157-290] days Of the 584 veins that were remapped, 72% (418/584) were durably isolated

Table 1. *Kaplan-Meier estimate with a median follow-up of 365 [323-386]


Conclusions

- Procedure times were short (58 min) despite a large number of operators with varied experience and workflow.
- There was a low rate of safety events (3.6%) and promising one-year efficacy rate (74%) in a large spectrum of AF patients.
- Similar one-year AF/AT recurrence rates among users with varying AF ablation experience demonstrated rapid adoption of the technology.
- Patients returning for a repeat ablation revealed a high rate of PVI with 72% of PVs durably isolated.



Integration with VersaCross will further simplify workflow

FAST TRACK to your therapy delivery in a **SINGLE SOLUTION**





VersaCross Brief Summary



VersaCross™ RF Wire

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross™ RF Wire is indicated for creation of an atrial septal defect in the heart.

CONTRAINDICATIONS: The VersaCross™ RF Wire is not recommended for use with any conditions that do not require the creation of an atrial septal defect. The Connector Cable is not recommended for use with any other Baylis RF Generator or any other device.

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during RF puncture procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. • The VersaCross™ RF Wire and Connector Cable are intended for single patient use only. Do not attempt to sterilize and reuse either devices. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another. Reuse may result in patient complications. • The VersaCross™ RF Wire must be used with the Connector Cable provided. Attempts to use it with other connector cables can result in electrocution of the patient and/or operator. • Do not use the VersaCross™ RF Wire with electrocautery or electrosurgery generators, connector cables or accessories as attempted use can result in patient and/or operator injury. • The Connector Cable must only be used with the RFP-100A Baylis RF Generator and the included VersaCross™ RF Wire. Attempts to use it with other RF Generators and devices can result in electrocution of the patient and/or operator. • The VersaCross™ RF Wire must be used with 0.035" compatible transseptal sheath and/or dilator devices. Use of incompatible accessory devices may damage the integrity of the VersaCross™ RF Wire or accessory devices and may cause patient injury. • The VersaCross™ RF Wire has only been validated for transseptal puncture use through VersaCross™ dilators which have been demonstrated to provide the required support for optimal function. • The VersaCross™ RF Wire is not intended for use with neonatal patients (i.e. less than one month of age). Do not attempt to treat neonatal patients with the VersaCross™ RF Wire.

PRECAUTIONS: • In order to prevent the risk of ignition, ensure that flammable materials are not present in the room during RF power application. • Careful manipulation of the VersaCross™ RF Wire must be performed to avoid vessel trauma. If resistance is encountered, DO NOT use excessive force to advance or withdraw the VersaCross™ RF Wire or ancillary sheath and/or dilator assembly. Excessive force may lead to bending or kinking of the device limiting advancement and retraction of sheath and/or dilator device. • The Baylis RF Generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the VersaCross™ RF Wire and/or DIP electrode, particularly when operating the device. • During power delivery, the patient should not be allowed to come in contact with ground metal surfaces. • If using electroanatomical mapping guidance, it is recommended to use it along with alternative imaging modality in the event there is loss of visibility of the device. • During power delivery, the patient should not be allowed to come in contact with ground metal surfaces.

ADVERSE EVENTS: Adverse events that may occur while creating an atrial septal defect include: • Tamponade • Sepsis/Infection • Thromboembolic episodes • Vessel perforation • Atrial Fibrillation • Myocardial Infarction • Vessel spasm • Sustained arrhythmias • Atrial Flutter • Hemorrhage • Vascular thrombosis • Perforation of the myocardium • Hematoma • Allergic reaction to contrast medium • Ventricular Tachycardia •

Pain and Tenderness • Arteriovenous fistula • Pericardial effusion • Tachycardia • Vascular Trauma • Additional Surgical Procedure • Wire entrapment/entanglement • Foreign body/wire fracture

EP-1504711-AA

VersaCross Connect™ Transseptal Dilator

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross Connect™ Transseptal Dilator is indicated for use in procedures where access to the left atrium via the transseptal technique is desired. **CONTRAINDICATIONS:** There are no known contraindications for this device.

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended. • The VersaCross Connect™ Transseptal Dilator and accompanying guidewire are intended for single patient use only. Do not attempt to sterilize and reuse the VersaCross Connect™ Transseptal Dilator or accompanying guidewire. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to do so may result in patient complications. • Careful manipulation must be performed to avoid cardiac damage or tamponade. Dilator advancement should be performed under imaging guidance. Fluoroscopic or echocardiographic imaging is recommended. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device. • DO NOT attempt to insert or retract the guidewire through a metal cannula or a percutaneous needle, which may damage the guidewire and may cause patient injury.

PRECAUTIONS: • The VersaCross Connect™ Transseptal Dilator is compatible with introducer sheaths 12.5Fr or larger. • The VersaCross Connect™ Transseptal Dilator is for use with specified models of 12F ID WATCHMANTM Access Sheath that are 75cm in length. • The VersaCross Connect™ Transseptal Dilator is compatible with 0.035" transseptal devices and guidewires or smaller. • The VersaCross Connect™ Transseptal Dilator is NOT compatible with transseptal needles such as the "NRG™ Transseptal Needle".

ADVERSE EVENTS: Adverse events that may occur while using the VersaCross Connect™ Transseptal Dilator and accompanying guidewire include: • Infection • Air embolus • Local nerve damage • Vessel trauma • Vessel spasm • Pseudoaneurysm • AV fistula formation • Atrial septal defect • Arrhythmias • Perforation and/or tamponade • Hematoma • Hemorrhage • Catheter entrapment • Embolic events • Valve damage •

Pericardial/pleural effusion

EP-1506005-AA

VersaCross™ Large Access Transseptal Dilator

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross™ Large Access Transseptal Dilator is indicated for use in procedures where access to the left atrium via the transseptal technique is desired. **CONTRAINDICATIONS:** There are no known contraindications for this device.

WARNINGS: Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended. • The VersaCross™ Large Access Transseptal Dilator and accompanying guidewire are intended for single patient use only. Do not attempt to sterilize and reuse the VersaCross™ Large Access Transseptal Dilator or accompanying guidewire. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to do so may result in patient complications. • Care should be taken to ensure that all air is removed from the dilator before infusing through the proximal hub. • Do not attempt direct percutaneous insertion of the dilator without a guidewire as this may cause vessel injury.

PRECAUTIONS: • Careful manipulation must be performed to avoid cardiac damage, or tamponade. Dilator and guidewire advancement should be performed under imaging guidance, such as fluoroscopy or echocardiography. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device. • The VersaCross™ Large Access Transseptal Dilator is compatible with introducer sheaths 12.5Fr or larger. • The VersaCross™ Large Access Transseptal Dilator is compatible with .035" transseptal devices and guidewires. • The VersaCross™ Large Access Transseptal Dilator is NOT compatible with transseptal needles such as the "NRG™ Transseptal Needle". • Do not reshape distal tip or curve of the guidewire. Excessive bending or kinking of the distal curve may damage the integrity of the wire or coating and lead to patient injury.

ADVERSE EVENTS: Adverse events that may occur while using the VersaCross™ Large Access Transseptal Dilator and accompanying guidewire include: • Infection • Air embolus • Local nerve damage • Vessel trauma • Vessel spasm • Pseudoaneurysm • AV fistula formation • Atrial septal defect • Arrhythmias • Perforation and/or tamponade • Hematoma • Hemorrhage • Catheter entrapment • Embolic events • Valve damage •

Pericardial/pleural effusion

EP-1506604-AA

VersaCross™ Steerable Sheath

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross™ Steerable Sheath kit is indicated for introducing various cardiovascular catheters to the heart, including the left side of the heart through the interatrial septum. **CONTRAINDICATIONS:** There are no known contraindications for this device.

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended. • The VersaCross™ Steerable Sheath kit is intended for single patient use only. Do not attempt to sterilize and reuse the VersaCross™ Steerable Sheath kit. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. • Do not attempt direct percutaneous insertion of the sheath without the dilator as this may cause vessel injury. • Maintain continuous hemodynamic monitoring throughout procedure • Provide continuous heparinized saline infusion while the introducer remains in vessel. • DO NOT attempt to insert or retract the guidewire through a metal cannula or a percutaneous needle, which may damage the guidewire and may cause patient injury.

PRECAUTIONS: • Careful manipulation must be performed to avoid cardiac damage, or tamponade. Sheath, dilator and guidewire advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device. • Avoid deflecting distal end of sheath during delivery and removal, otherwise damage to vessels may occur. • The VersaCross™ Steerable Sheath kit is not compatible with transseptal needles such as the "NRG™ Transseptal Needle".

ADVERSE EVENTS: Adverse events that may occur while using the VersaCross™ Steerable Sheath include: • Infection • Air embolus • Local nerve damage • Vasovagal reaction • Dissection • Vessel spasm • AV fistula formation • Atrial septal defect • Pseudoaneurysm • Aortic puncture • Arrhythmias • Perforation and/or tamponade • Hematoma • Hemorrhage • Catheter entrapment • Embolic events • Stroke • Valve damage •

Myocardial Infarction • Pericardial/pleural effusion • Pulmonary edema • Coronary artery spasm and/or damage • Vessel trauma • Pacemaker/defibrillator lead displacement

EP-1506705-AA

VersaCross™ Transseptal Sheath

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross™ Transseptal Sheath kit is used for the percutaneous introduction of various types of cardiovascular catheters and guidewires to all heart chambers, including the left atrium via transseptal perforation/puncture. **CONTRAINDICATIONS:** There are no known contraindications for this device.

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. The use of echocardiography is recommended. • The VersaCross™ Transseptal Sheath kit is intended for single patient use only. Do not attempt to sterilize and reuse the VersaCross™ Transseptal Sheath kit. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another. • Do not attempt direct percutaneous insertion of the sheath without the dilator as this may cause vessel injury. • Careful manipulation must be performed to avoid cardiac damage or tamponade. Sheath advancement should be done under fluoroscopic guidance. Echocardiographic guidance is also recommended.

PRECAUTIONS: • Careful manipulation must be performed to avoid cardiac damage, or tamponade. Sheath, dilator and guidewire advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device. • The VersaCross™ Transseptal Sheath is compatible with introducer sheaths 11Fr or larger. • The VersaCross™ Transseptal Sheath and Dilator are compatible with transseptal devices and guidewires .035" or smaller. • The VersaCross™ Transseptal Sheath kit is NOT compatible with transseptal needles such as the "NRG™ Transseptal Needle".

ADVERSE EVENTS: Adverse events that may occur while using the VersaCross™ Transseptal Sheath kit include: • Infection • Air embolus • Local nerve damage • Hemorrhage • Embolic events • Vessel spasm • AV fistula formation • Atrial septal defect • Pseudoaneurysm • Perforation and/or tamponade • Arrhythmias • Pericardial/pleural effusion • Hematoma • Vessel trauma • Valve damage • Catheter entrapment

EP-1506605-AA

VersaCross™ Transseptal Dilator

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Instructions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INDICATIONS FOR USE: The VersaCross™ Transseptal Dilator is indicated for introducing various cardiovascular catheters to the heart, including the left side of the heart through the interatrial septum. **CONTRAINDICATIONS:** There are no known contraindications for this device.

WARNINGS: • Laboratory staff and patients can undergo significant x-ray exposure during interventional procedures due to the continuous usage of fluoroscopic imaging. This exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects. Therefore, adequate measures must be taken to minimize this exposure. • The VersaCross™ Steerable Sheath kit is intended for single patient use only. Do not attempt to sterilize and reuse the VersaCross™ Steerable Sheath kit. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to follow this instruction may result in patient complications • Maintain continuous hemodynamic monitoring throughout procedure • Provide continuous heparinized saline infusion while the introducer remains in vessel.

PRECAUTIONS: • Careful manipulation must be performed to avoid cardiac damage, or tamponade. Sheath, dilator and guidewire advancement should be done under fluoroscopic guidance. If resistance is encountered, DO NOT use excessive force to advance or withdraw the device. • The VersaCross™ Steerable Sheath kit is intended for single patient use only. Do not attempt to sterilize and reuse the VersaCross™ Steerable Sheath kit. Reuse can cause the patient injury and/or the communication of infectious disease(s) from one patient to another. Failure to follow this instruction may result in patient complications • Maintain continuous hemodynamic monitoring throughout procedure • Provide continuous heparinized saline infusion while the introducer remains in vessel.

ADVERSE EVENTS: Adverse events that may occur while using the VersaCross™ Sheath include: • Infection • Air embolus • Local nerve damage • Vasovagal reaction • Dissection • Vessel spasm • AV fistula formation • Atrial septal defect • Pseudoaneurysm • Aortic puncture • Arrhythmias • Perforation and/or tamponade • Hematoma • Hemorrhage • Catheter entrapment • Thromboembolic events • Stroke • Valve damage •

Myocardial Infarction • Pacemaker/defibrillator lead displacement • Pulmonary edema • Coronary artery spasm and/or damage • Vessel trauma • Pericardial/pleural effusion

EP-1506213-AA