



Go Beyond Surveillance Biopsies with HeartCare

The first and only multi-modal test with studies demonstrating a significant reduction in surveillance biopsies and revolutionizing post heart transplant care



Sam D., heart transplant recipient, and his wife

Unmet Need

Heart Transplant Surveillance is Complex

Risk of long-term complications require a holistic approach

TOO MUCH IMMUNOSUPPRESSION


Increased risk of:

- Infections
- Kidney injury
- Malignancy
- Post-transplant diabetes

TOO LITTLE IMMUNOSUPPRESSION

Increased risk of:


- Rejection leading to potential failure of the organ
- Cardiac Allograft Vasculopathy (CAV)



Unmet Need

Traditional Tools for Surveillance Monitoring May Not Allow for Early Intervention

This may impact long term outcomes




Unmet Need

Traditional Biopsies Have Significant Limitations for Rejection Surveillance

- Resource intensive
- Risk of exposure to pathogens
- Sampling errors and limitations
- Invasive, painful, and frequent (sometimes >10 Biopsies in the 1st year post-transplant)
- Interobserver variability in interpretation – for example:

Almost 40% of rejections called by local pathologists were downgraded when reviewed by a panel of central pathologists

Agreement between local and central pathology for rejection ¹		
	Central Grades	
Local Grades	2R Rejection	< 2R Rejection
2R Rejection	46/76 (60.5%)	30/76 (39.5%)
< 2R	63/862 (7.3%)	799/862 (92.7%)



Introducing HeartCare

HeartCare is a non-invasive leading indicator of graft injury and immune activation/quiescence

Unmet Need

Introducing HeartCare


HeartCare as a Surveillance Tool

Monitoring Rejection and Beyond with HeartCare

Clinical Studies

Logistics Data

Your Partner in Transplant Care



Introducing HeartCare


HeartCare Is a Multi-Modality Approach that Utilizes Two Complementary Technologies



Immune Quiescence	Detection	
	ACR	AMR
+++	+++	+++



Gene Expression Profiling



Genes in the AlloMap Signature

- Lymphocyte Activation: **SEM7A**
- Cell Migration: **RHOU**
- T Cell Priming: **PDCD1, ITGA4**
- Inflammation (Hematopoietic Proliferation): **MARCH8, WDR40A**
- Steroid Sensitive: **IL1R2, FLT3, ITGAM**
- Platelet Activation: **PF4, C6orf25**


Immune Quiescence	Detection
	ACR
+++	++

+



Donor-Derived Cell-Free DNA

dd-cfDNA No Graft Injury **dd-cfDNA During Graft Injury**



Detection	
ACR	AMR
++	+++

Introducing HeartCare

AlloMap Heart is an Established Standard in Heart Transplant

Unmet Need

Introducing HeartCare

HeartCare as a Surveillance Tool

Monitoring Rejection and Beyond with HeartCare

Clinical Studies

Logistics Data

Your Partner in Transplant Care



Used in 90% of US heart transplant centers and 1 in 2 new patients



Guideline recommended in 2010 and 2023*



Extensively validated including publication in NEJM



COVERED BY MEDICARE**

**When coverage criteria is met



FDA CLEARED IN 2008 & 2023

FDA 510(k)-cleared to identify recipients at low probability of ACR ≥55 days post-transplant and ≥15 years old

* References to ISHLT are offered solely to support AlloMap's FDA indications and should not be construed as supporting any other use.

Introducing HeartCare

HeartCare Provides Peace-of-Mind During Surveillance

Unmet Need

Introducing HeartCare


HeartCare as a Surveillance Tool

Monitoring Rejection and Beyond with HeartCare

Clinical Studies

Logistics Data

Your Partner in Transplant Care



[^] Post-transplant Period $\geq 2-6$ months
^{^^} Post-transplant period ≥ 6 months

* clinical correlation is required
** as defined by their institutional protocol

Introducing HeartCare

New ISHLT Guidelines² Support Use of HeartCare Solutions, AlloMap Heart and AlloSure Heart, in Routine Monitoring of Heart Transplant Patients

The new 2023 ISHLT Guidelines mention the following:

2010 Guideline Recommendation	2023 Update Guideline Recommendation
Earlier use of AlloMap Heart starting at two months post-transplant based on strength of clinical studies.	
Gene Expression Profiling (Allomap) can be used to rule out the presence of ACR of grade 2R or greater in appropriate low-risk patients, starting at 6 months after HT. <i>Class IIa, Level of Evidence: B</i>	Gene Expression Profiling (GEP) (i.e., Allomap) of peripheral blood can be used in low-risk patients starting at 2 months after HT to identify adult recipients who have low risk of current ACR to reduce the frequency of EMB. <i>Class IIa, Level of Evidence: B</i>
AlloMap gene expression profiling (GEP) and donor-derived cell-free DNA (dd-cfDNA) for routine post-transplant monitoring.	
No mention of non-invasive monitoring in the follow up visits to monitor rejection	Follow-up visits to monitor for rejection may include non-invasive rejection monitoring - Gene Expression Profiling [AlloMap] and donor-derived cell-free DNA <i>Class I, Level of Evidence B.</i>
Remote use of GEP and dd-cfDNA heart transplant surveillance (HeartCare), as evidenced during the COVID pandemic.	
No recommendation	Patient Management During a Pandemic Efforts should be made to reduce visits by clinically stable heart transplant patients to medical facilities by shifting blood testing to the patients' homes. Remote drawing of blood samples can include screening tests to determine if patients require endomyocardial biopsies using gene expression profiling and donor derived cell-free DNA assays. <i>Class I, Level of Evidence C</i>

CareDx is the **ONLY** provider of **Gene Expression Profiling (AlloMap Heart)** and **Donor-Derived Cell-Free DNA (AlloSure Heart)** for heart transplant recipients

References to ISHLT are offered solely to support AlloMap's FDA indications and should not be construed as supporting any other use. AlloMap should be solely used in conjunction with standard clinical assessment

HeartCare as a Surveillance Tool

HeartCare Offers a Synergistic Approach to Post-transplant Surveillance



A Non-Invasive Surveillance Tool

- Enhanced ACR identification that is superior to using one test
- Significant Reduction in Surveillance Biopsies (confirmed by multiple studies)




Measures mRNA of the recipient



Measures immune response



Originates from WBCs



Measures dd-cfDNA from the transplanted heart




Originates from the graft

HeartCare as a Surveillance Tool

High HeartCare Results are Associated with a Higher Chance of ACR

Results from two independent studies^{3,4} demonstrate increased ACR with concordantly elevated HeartCare results




High AlloMap is: ≥ 30 for ≥ 2 -6 months or ≥ 34 for ≥ 6 months | High AlloSure is $\geq 0.2\%$

HeartCare as a Surveillance Tool

High HeartCare Results Increased the Odds of ACR by ~14X

Data derived from the combination of two independent single center studies^{3,4}

The chance of a biopsy revealing ACR is **greater with a high HeartCare** than with either a high AlloMap or a high AlloSure alone




Positive Likelihood Ratio gives the change in the odds of having a diagnosis in patients with a positive test.

A LR+ is mathematically defined as sensitivity / (1-specificity)

High AlloMap is: ≥ 30 for $\geq 2-6$ months or ≥ 34 for ≥ 6 months | High AlloSure is $\geq 0.2\%$

HeartCare as a Surveillance Tool

ACR Surveillance with HeartCare – Proposed Management^{1,3,4,5}



HeartCare as a Surveillance Tool

Clinical Interpretation for HeartCare - ACR Surveillance

HeartCare helps clinicians precisely identify patients who are at higher risk of ACR

<p>HIGH ALLOMAP / LOW ALLOSURE</p> <p>A biopsy is unlikely to reveal ACR</p> <p>Consider biopsy or repeat HeartCare testing earlier if:</p> <ul style="list-style-type: none"> • AlloMap elevation is new and AlloSure level has risen or is near threshold <p>Consider other pathological causes of an increased AlloMap:</p> <ul style="list-style-type: none"> • CMV infection <p style="text-align: right; font-size: small;">Medication adherence review recommended</p>	<p>HIGH ALLOMAP / HIGH ALLOSURE</p> <p>Relatively high probability that ACR is present</p> <ul style="list-style-type: none"> • Biopsy should be considered <p style="text-align: right; font-size: small;">Medication adherence review recommended</p>
<p>LOW ALLOMAP / LOW ALLOSURE</p> <p>Low Likelihood of ACR</p> <ul style="list-style-type: none"> • A biopsy is unlikely to reveal ACR • Continue protocol immuno-optimization 	<p>LOW ALLOMAP / HIGH ALLOSURE</p> <p>A biopsy is unlikely to reveal ACR</p> <p>Consider a biopsy or repeat HeartCare testing earlier if:</p> <ul style="list-style-type: none"> • Recent treatment for rejection (<21 days) or current prednisone >20mg • At risk of Antibody Mediated Rejection / markedly elevated AlloSure level • AlloSure elevation is new and AlloMap level has risen or is near threshold <p>Consider other possible pathological causes of an increased AlloSure:</p> <ul style="list-style-type: none"> • Cardiac allograft vasculopathy • Severe infection • Antibody Mediated Rejection (AMR) / Donor specific antibodies <p style="text-align: right; font-size: small;">Medication adherence review recommended</p>



The table is provided for informational purposes only and is not intended as medical advice. A physician's test selection and interpretation, diagnosis, and patient management decisions should be based on his/her education, clinical expertise, current guidelines, and assessment of the patient.

This table is designed for the context of surveillance testing for ACR; For patients that are at risk of AMR or being tested in other clinical context, different guidance may apply.

High AlloMap is ≥ 30 for $\geq 2-6$ months or ≥ 34 for ≥ 6 months | High AlloSure is $\geq 0.2\%$

HeartCare as a Surveillance Tool

Two Independent Studies^{4,6} Demonstrate Significant Reduction In Surveillance Biopsies When Using HeartCare




HeartCare as a Surveillance Tool

HeartCare Allows for Reduction in Biopsies Without Impacting the Outcomes vs. GEP Alone


This independent study⁶ compares paired use of AlloSure dd-cfDNA and AlloMap GEP testing against GEP testing alone

The results were similar to the post-transplant survival, rejection free survival and graft function at 1 year, while requiring significantly fewer endomyocardial biopsies

Total Biopsies in the First Year Post Transplant



First Year Rejection Free Survival from Acute Cellular Rejection 2R/3R by Protocol



HeartCare as a Surveillance Tool

Clinicians Rely on HeartCare When Making Clinical Decisions

An independent study⁶ shows that clinicians find individual value in AlloMap and AlloSure, and use both tests for patient management

**In patients with a High AlloSure and a Low AlloMap,
More than 75% of ACR Surveillance biopsies were deferred**

Only AlloSure Heart

If 10 Patients had a High AlloSure



10 BIOPSIES

HeartCare Testing

If 10 Patients had a High AlloSure with a Low AlloMap



~2 BIOPSIES

~8 Biopsies Deferred


Data from Henricksen et al⁶

AlloSure Heart Has Been Clinically Validated to Differentiate Rejection vs. No Rejection

AlloSure Heart was clinically validated in the multi-center prospective D-OAR study⁵

Clinical Validation:

- **dd-cfDNA levels are over 2X greater in the rejection groups vs. the no rejection group**
- **No rejection vs. all rejection, p= 0.005**





- ACR = Grade 2R and Grade 3R
- Grades 0R and 1R are in the No Rejection group
- AMR = pAMR1 and pAMR2

- Unmet Need
-
- Introducing HeartCare
-
- HeartCare as a Surveillance Tool
-
- Monitoring Rejection and Beyond with HeartCare
-
- Clinical Studies
-
- Logistics Data
-
- Your Partner in Transplant Care

Monitoring Rejection and Beyond with HeartCare


AlloMap's Surveillance was Compared to Biopsy in the Landmark NEJM Paper⁷

This was followed by the ISHLT-cited eIMAGE randomized trial⁸ in which AlloMap was compared to biopsy beginning at **55 days** after transplant.



No. at Risk	0	100	200	300	400	500	600	700	800
Routine biopsies	305	278	252	221	181	160	137	137	73
Gene-expression profiling	297	273	252	207	177	162	133	130	36

Kaplan–Meier Estimates of the Time to the Composite Primary Outcome (in IMAGE):
 This image shows the time to the first occurrence of any of the following primary events: rejection with hemodynamic compromise, graft dysfunction due to other causes, death, or retransplantation. Only the first event that was part of the composite primary outcome was considered.




eIMAGE: Kaplan–Meier curve demonstrating probability of eventfree progression for GEP and EMB groups in the first 18 months post transplant

Monitoring Rejection and Beyond with HeartCare


AlloSure Heart Can Provide Information on DSA and CAV

REJECTION


Patients with AlloSure Heart levels >0.15% were more likely to develop De Novo DSAs within year 1 post-transplant⁹



Elevated AlloSure Levels are an Independent Predictor for dnDSA (p<0.001)¹⁰



Elevated AlloSure Levels are Associated with Cardiac Allograft Vasculopathy¹¹



Significant CAV was present in 63% of patients with high levels of dd-cfDNA compared to 35% in the low dd-cfDNA group (p=0.047)

Unmet Need

Introducing HeartCare

HeartCare as a Surveillance Tool

Monitoring Rejection and Beyond with HeartCare

Clinical Studies

Logistics Data

Your Partner in Transplant Care

Clinical Studies

HeartCare Has Been Clinically Validated in Multi-Center Prospective Studies Including >4,500 Patients

2006	Cardiac Allograft Rejection Gene Expression Observation (CARGO) First multicenter validation of GEP test to identify heart txp at low risk mod/severe rejection	<i>AJT</i>	N = 222
2010	Invasive Monitoring Attenuation through Gene Expression (IMAGE) Prospective, randomized, controlled, multi-center trial where Surveillance for rejection with AlloMap (gene expression profiling) was compared to routine biopsy (standard of care)	<i>NEJM</i>	N = 602
2015	Early Invasive Monitoring Attenuation through Gene Expression (eIMAGE) Randomized controlled trial comparing GEP to EMB as early as 55 days post transplant	<i>Circ: Heart Fail.</i>	N = 60
2016	CARGO II The CARGO II study validated the CARGO I study findings of GEP score performance in a predominantly European based cohort	<i>Eur Heart J.</i>	N = 472
2019	Outcomes AlloMap® Registry (OAR) Registry 5-year follow-up of patients from 2013-2019	<i>JHLT</i>	N = 1,504
2019	Donor-derived cell free DNA Outcomes AlloMap Registry (D-OAR) Subset of patients co-enrolled in OAR Prospective, Multi-center Registry with 26 Centers to determine whether dd-cfDNA level in a HT recipients' blood can differentiate rejection from no rejection	<i>AJT</i>	N = 740
2021	Surveillance HeartCare Outcomes Registry (SHORE) Prospective registry aimed to capture events in patients during their first 5 years post-transplant	In progress	N = 2707



Clinical Studies

The Evolution of AlloMap → AlloSure → HeartCare



Invasive Monitoring Attenuation through Gene Expression (IMAGE) was a prospective, randomized, controlled, multi-center trial⁷

Surveillance for rejection with AlloMap® (gene expression profiling) was compared to routine surveillance biopsy (standard of care)



- Unmet Need
- Introducing HeartCare
- HeartCare as a Surveillance Tool
- Monitoring Rejection and Beyond with HeartCare
- Clinical Studies
- Logistics Data
- Your Partner in Transplant Care

Clinical Studies

The Evolution of AlloMap → AlloSure → HeartCare




Early Invasive Monitoring Attenuation through Gene Expression (eIMAGE) was a single center, randomized, parallel 2-arm study beginning at **55 days** post heart transplant.⁸

- The primary endpoints between GEP and EMB groups were compared (p= .45)
- Components of the composite primary endpoints were: first year IVUS results and 1 year echo LV ejection fraction
- Patient satisfaction score for method of surveillance at 1-year was higher in the GEP vs. EMB group (p = .003)

Primary endpoint (composite): death or re-transplantation, hemodynamic compromise with rejection, graft dysfunction (hemodynamic compromise without rejection)

Co-primary endpoint: first-year IVUS changes post-transplant



Kaplan–Meier curve demonstrating probability of eventfree progression for GEP and EMB groups in the first 18 months post transplant


Clinical Studies

The Evolution of AlloMap → AlloSure → HeartCare



Donor-Derived Cell-Free DNA-**Outcomes** AlloMap Registry (**D-OAR**) was a prospective, multi-center registry with 26 centers to determine whether dd-cfDNA level in a HT recipients' blood can differentiate rejection from no rejection.⁵

AlloSure levels were over 2x greater in the rejection groups vs. the no rejection group



ACR = 2R and 3R; AMR = pAMR1 and pAMR2

HeartCare - Summary

Proven multi-modality surveillance for heart transplant recipients.

**HeartCare can be used to provide
a comprehensive overview of graft injury
and immune quiescence.**



HeartCare can be used for post-transplant surveillance to identify stable patients who are at a higher than anticipated risk of current rejection



HeartCare can also better distinguish between the subset of the population with a higher probability of rejection than either AlloMap or AlloSure Heart alone



Two independent publications have demonstrated that with HeartCare, the avoidance of biopsies is superior to either AlloMap or AlloSure Heart alone


Logistics Data

HeartCare Can Be Used For Longitudinal Surveillance

The Heart Allograft Routine Testing Schedule (HARTS)

The Heart Allograft Routine Schedule is recommended based upon the use of AlloMap in the IMAGE and eIMAGE trials and use in clinical practice. All testing should be performed when medically necessary, in accordance with a physician's guidance. The SHORE registry study is assessing the clinical utility of this testing interval.

HARTS Timeline



*Under FDA-cleared labeling, AlloMap is indicated for use in heart transplant patients who are ≥15 years old and ≥55 days post transplant. All specimens submitted for testing and billed to third party payers must support the medical necessity of testing when a test is ordered consistent with clinical judgment and applicable payer policies <https://clinicaltrials.gov/ct2/show/NCT03695601>

Real-World HeartCare Surveillance*

The Stanford University heart transplant protocol for rejection surveillance

Time	Surveillance	Donor Specific Antibody	Echocardiogram	Right Heart Catheterization
Week 2	Biopsy	If antibody crossed	—	Yes
Week 3	Biopsy	If antibody crossed	—	Yes
1 month	AlloSure/Biopsy	Yes	Limited	Yes
Week 6	AlloSure/Biopsy	—	—	Yes
2 months	HeartCare or Biopsy	—	Limited	Yes
3 months	HeartCare or Biopsy	Yes	—	Yes
4 months	HeartCare or Biopsy	—	Limited	—
5 months	HeartCare or Biopsy	—	—	—
6 months	HeartCare or Biopsy	Yes	Limited	—
8 months	HeartCare or Biopsy	—	—	—
10 months	HeartCare or Biopsy	Yes	—	—
12 months	HeartCare or Biopsy	Yes	Full or DSE if no angiogram	Yes
15 months	HeartCare or Biopsy	—	—	—
18 months	HeartCare or Biopsy	—	Limited	—
21 months	HeartCare or Biopsy	—	—	—
24 months	HeartCare or Biopsy	Yes	Full or DSE if no angiogram	Yes (if angiogram)

*This protocol is provided as an example only. All clinical decisions should be made directly by the referring provider.
DSE: Dobutamine Stress Echocardiogram

- Unmet Need
-
- Introducing HeartCare
-
- HeartCare as a Surveillance Tool
-
- Monitoring Rejection and Beyond with HeartCare
-
- Clinical Studies
-
- Logistics Data
-
- Your Partner in Transplant Care

Logistics Data

Sample Report - HeartCare

- Individual test results (AlloMap or AlloSure Heart) are released when available on individual reports
- Combined HeartCare results are provided on an integrated report that includes both AlloMap and AlloSure Heart when both test results are available

Patient Name: _____

Unique Patient Identifier: _____

Date of Birth: _____

Transplant Date: _____

Client: _____

Ordered By: _____

ALLOMAP[®] RESULTS

Accession ID: _____

Report ID: _____

Client Specimen ID: _____

Specimen Date: _____

Date Received: _____

Date Reported: _____

AlloMap Score

33

AlloMap Score Variability

< 0.04%

Improved reportable range 0.04% - 16% (previously 0.12% - 16%).

Thresholds associated with rejection have NOT changed.

Date	Score
18 Jul 2023	25
22 Aug 2023	31
11 Sep 2023	33
12 Oct 2023	24
09 Nov 2023	26
07 Dec 2023	33

NR - No Result, US - Unacceptable Sample
Refer to the individual AlloMap and AlloMap Variability Reports when using these results as part of overall patient management.

ALLOSURE HEART[®] RESULTS

Accession ID: _____

Report ID: _____

Client Specimen ID: _____

Specimen Date: _____

Date Received: _____

Date Reported: _____

AlloSure Heart Result

< 0.04%

Improved reportable range 0.04% - 16% (previously 0.12% - 16%).

Thresholds associated with rejection have NOT changed.

Date	Result
18 Jul 2023	+0.04
22 Aug 2023	+0.04
11 Sep 2023	+0.04
12 Oct 2023	+0.04
09 Nov 2023	+0.04
07 Dec 2023	+0.04

INTERPRETATION OF ALLOSURE HEART RESULTS

AlloSure Heart measures the percent of donor-derived cell-free DNA (dd-cfDNA) in the total cell-free DNA present in the heart transplant recipients.¹ When interpreting AlloSure Heart:

- A result $\geq 0.15\%$ dd-cfDNA is associated with a higher probability of acute cellular rejection (ACR) and/or antibody mediated rejection (AMR).¹ See the HeartCare interpretation guide at caredx.com/heartcare for information on AlloSure Heart performance at other cutoffs.
- A result $< 0.15\%$ dd-cfDNA is associated with a lower risk of ACR and/or AMR with a $> 97\%$ negative predictive value.¹ 76% of plasma samples from heart transplant patients without clinical signs or symptoms of allograft dysfunction were observed to have a dd-cfDNA value $< 0.15\%$.

Clinical validity of AlloSure Heart was established in heart transplant recipients age 15 or older and at least 14 days post-transplant. AlloSure Heart should be interpreted in the context of clinical findings and relevant patient history. AlloSure Heart is not intended for recipients of multi-organ transplants, recipients who are pregnant, or recipients of allogeneic bone marrow transplant. Transfusion of blood components containing white blood cells in the 30 days prior to blood draw may lead to an aberrant result. Endomyocardial biopsies performed in the 24 hours prior to AlloSure Heart specimen collection may elevate dd-cfDNA. For more information about AlloSure Heart, please visit CareDx.com/heartcare.

¹ Grskovic et al., J Mol Diagn 2016; ² Khush et al., Am J Transplant 2019.

The AlloSure Heart Donor-Derived Cell-Free DNA Test was developed and its performance characteristics were determined by the CareDx laboratory. The test has not been cleared or approved by the U.S. Food and Drug Administration; FDA approval is not required. The laboratory is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA '88) and accredited by the College of American Pathologists (CAP) as qualified to perform high complexity clinical laboratory testing.

The contents of this report are confidential and intended solely for the use of authorized personnel.

3260 Bayshore Blvd., Brisbane, California 94005 USA
P: 888-255-6627 F: 415-287-2456
www.caredx.com | www.alloasure.com

Lab Director:
William Patrick Joseph, M.D.
CLIA #: 05D1029609 | CAP #: 7193868

F-T-CF-0012 Revision 7 Effective 2023-04

28

Ordering HeartCare

HeartCare may be ordered in several ways:



CarePortal

HeartCare may be ordered electronically, via the CarePortal. One-time tests or custom orders are available through the portal.



HeartCare Test Requisition Form (TRF)

A paper test requisition form is available for ordering HeartCare. Please contact Customer Care for delivery of the hardcopy HeartCare TRFs



Direct EMR Integration

HeartCare may be ordered directly from a hospital EMR by integrating with CareDx's lab software

Your Partner in Transplant Care

Solutions Throughout the Transplant Journey

A suite of innovative transplant management solutions, enabling more cohesive care



Operational Support and Quality Improvement



Referral Management

Connects referring provider, transplant center, and patient data



Customizable Administrative Service

Manages specific tasks for a transplant center



Quality Management Software

Aggregates data and leverages predictive analytics to help transplant centers improve quality metrics



Informed Decision Support



Informed Decisions

A suite of AI-derived offerings that leverage multiple clinical inputs to assist clinicians in making more informed decisions



Donor-Derived, Cell-Free DNA

Provides early signal of allograft injury



AlloMap + AlloSure

Provides insights into allograft injury and immune quiescence via dd-cfDNA and gene expression profiling



Patient Management



Specialty Pharmacy

Provides comprehensive pharmacy support via Patient Care Teams to increase adherence



Medication Management and Education

Simplifies complex medication regimens to increase adherence



Remote Patient Monitoring

Provides real-time, at home vital sign measurements

Your Partner in Transplant Care

CareDx is Your Dedicated Partner in Transplant Care

CareDx supports our providers and patients by providing dedicated Patient Care Managers to facilitate appointment scheduling and blood draw logistics at no additional cost.*

Services

Patient Care Managers


Our Patient Care Managers include nurses and former transplant coordinators who understand the complexity in managing heart transplant patients.

RemoTraC®

CareDx provides mobile phlebotomy services for patient tests. There are no additional out of pocket costs beyond any coinsurance or deductible required by the patient's health plan.

MedActionPlan

CareDx offers a robust medication adherence platform to support transplant patients with better medication adherence to optimize their transplant outcomes



*CareDx provides specimen collection through mobile phlebotomists as part of its laboratory service so patients should not have additional out of pocket expenses beyond any coinsurance/deductibles required by their health plans for the test

References

1. Crespo-Leiro MG, Stypmann J, Schulz U, et al. Clinical usefulness of gene-expression profile to rule out acute rejection after heart transplantation: CARGO II. *Eur Heart J*. 2016;37(33):2591-2601. doi:10.1093/eurheartj/ehv682.
2. Velleca A, Shullo MA, Dhital K, et al. The International Society for Heart and Lung Transplantation (ISHLT) guidelines for the care of heart transplant recipients. *J Heart Lung Transplant*. 2023;42(5):e1-e141. doi:10.1016/j.healun.2022.10.015
3. Data Adapted from Rodgers N, Gerding B, Cusi V, et al. Comparison of two donor-derived cell-free DNA tests and a blood gene-expression profile test in heart transplantation. *Clin Transplant*. 2023;37(6):e14984. doi:10.1111/ctr.14984
4. Gondi KT, Kao A, Linard J, et al. Single-center utilization of donor-derived cell-free DNA testing in the management of heart transplant patients. *Clin Transplant*. 2021;35(5):e14258. doi:10.1111/ctr.14258
5. Khush KK, Patel J, Pinney S, et al. Noninvasive detection of graft injury after heart transplant using donor-derived cell-free DNA: A prospective multicenter study. *Am J Transplant*. 2019;19(10):2889-2899. doi:10.1111/ajt.15339. ; D-OAR is a sub-study of the Outcomes AlloMap Registry (OAR).
6. Henriksen EJ, Moayed Y, Purewal S, et al. Combining donor derived cell free DNA and gene expression profiling for non-invasive surveillance after heart transplantation [published online ahead of print, 2022 May 12]. *Clin Transplant*. 2022:e14699. doi:10.1111/ctr.14699.
7. Pham MX, Teuteberg JJ, Kfoury AG, et al. Gene-expression profiling for rejection surveillance after cardiac transplantation. *N Engl J Med*. 2010;362(20):1890-1900. doi:10.1056/NEJMoa0912965.
8. Kobashigawa J, Patel J, Azarbal B, et al. Randomized pilot trial of gene expression profiling versus heart biopsy in the first year after heart transplant: early invasive monitoring attenuation through gene expression trial. *Circ Heart Fail*. 2015;8(3):557-564. doi:10.1161/CIRCHEARTFAILURE.114.001658.
9. Teuteberg, J., Kobashigawa, J., Shah, P., Ghosh, S., Ross, D., DePasquale, E., & Khush, K. (2021). Donor-Derived Cell-Free DNA Predicts De Novo DSA after Heart Transplantation. *J Heart Lung Transplant*, 40(4), S30.
10. DePasquale E, Kobashigawa J, Pinney S, Khush K. dd-cfDNA as a Risk Factor for Initiating De-Novo Donor Specific Antibodies in Heart Transplantation [abstract]. *Am J Transplant*. 2020; 20 (suppl 3). <https://atcmeetingabstracts.com/abstract/dd-cfdna-as-a-risk-factor-for-initiating-de-novo-donor-specific-antibodies-in-heart-transplantation/>. Accessed August 3, 2022.
11. Holzhauser L, Clerkin KJ, Fujino T, et al. Donor-derived cell-free DNA is associated with cardiac allograft vasculopathy. *Clin Transplant*. 2021;35(3):e14206. doi:10.1111/ctr.14206.