

CONFIDENTLY TREAT WIDE NECKS WITH ESAR

EndoSuture Aneurysm Repair (ESAR) offers a durable, reinforced seal and protects against neck dilatation to minimize Type Ia endoleaks


What are the challenges after EVAR in wide necks?

Neck Dilatation

Neck diameter is an independent risk factor for neck dilatation.¹ Diameter $\geq 25\text{mm}$ has greater risk of dilatation after Endovascular Aneurysm Repair (EVAR) ($p=0.02$).²

Type Ia Endoleaks

Wide necks are more likely to develop Type Ia endoleaks ($p=0.049$);³ which transmit the highest pressures into the sac.⁴



Neck diameter predicts Type Ia endoleak risk*:

*observed with multiple endografts

GREAT Registry (n = 2678)

$\geq 25\text{mm}$ predicts greater Type Ia risk ($p=0.007$)⁵ **Device Included:** Excluder™ graft


ANCHOR Registry (n = 221)

$\geq 26\text{mm}$ (at lowest renal) predicts development of Type Ia ($p=0.013$)⁶ **Devices Included:** Endurant™ system, Excluder™, Zenith™, Talent™, AneuRx™ grafts

™*Third party brands are trademarks of their respective owners.

ESAR enhances the proximal seal zone competency over standard EVAR

A Creates a durable, reinforced seal⁷⁻⁹



STRONGER ATTACHMENT^{8,9}
Over EVAR⁸ alone, via secure transmural wall fixation¹⁰




Image courtesy of Dr. Giovanni Pratesi

RADIAL FIXATION¹¹
95.9% implants with adventitial penetration^{7*}
Radial support just as with sutures¹¹

*ANCHOR Primary AAA, at intended location

B Protects against neck dilatation¹



Benchmark data on file at Medtronic

LOSS OF APPPOSITION WITH EVAR
Neck expansion beyond nominal graft diameter^{1,12}



Benchmark data on file at Medtronic

MAINTAIN APPPOSITION WITH ESAR
Stability of a surgical anastomosis⁸
EndoAnchor™ implants protect where they are deployed¹

C Minimizes Type Ia endoleaks




Image courtesy of Dr. Jean-Paul De Vries

MINIMAL TYPE IA IN HOSTILE NECKS

	ANCHOR Primary AAA ¹³ n = 712	EVAR Only ¹⁴ n = 199
Type Ia	3.4% (4/117)	9.5% (19/199)
Time at 4 years	88.6%	100%
Hostile*		

Results from separate trials; head to head data may differ * hostile necks per SVS criteria




Image courtesy of Dr. Apostolos Tassiopoulos

MINIMAL TYPE IA REINTERVENTIONS
97.2% FF reinterventions for Type Ia endoleak (through 4 years)⁷

"EndoAnchor™ implants are an opportunity for EVAR improvement"¹⁵

Dr. William Jordan Jr., M.D.
Professor and Chief of Vascular Surgery, Emory University

Medtronic

REFERENCES

- ¹ Tassiopoulos AK, Monastiriotis S, Jordan WD, Muhs BE, Ouriel K, De Vries JP. Predictors of early aortic neck dilatation after endovascular aneurysm repair with EndoAnchors. *J Vasc Surg.* 2017;66(1):45-52
- ² Cao P, Verzini F, Parlani G, et al. Predictive factors and clinical consequences of proximal aortic neck dilatation in 230 patients undergoing abdominal aorta aneurysm repair with self-expandable stent-grafts. 2003 June; 37(6):1200-5.
- ³ Oliveira NFG, Golcalves FB, Ultee K, et al. Patients with large neck diameter have a higher risk of type 1A endoleaks and aneurysm rupture after standard endovascular aneurysm repair. *J Vasc Surg.* 2019 March;69(3):783-791.
- ⁴ Li J, Tian X, Wang Z, et al. Influence of endoleak positions on the pressure shielding ability of stent-graft after endovascular aneurysm repair (EVAR) of abdominal aortic aneurysm (AAA). *Biomed Eng Online.* 2016;15(Suppl 2): 135.
- ⁵ Howard et al. 2018. Influence of Proximal Aortic Neck Diameter on Durability of Aneurysm Sealing and Overall Survival in Patients Undergoing Endovascular Aneurysm Repair. Real World Data from the Gore Global Registry for Endovascular Aortic Treatment (GREAT)
- ⁶ Jordan et al. 2015. Outcome-based anatomic criteria for defining the hostile aortic neck
- ⁷ ANCHOR 4-yr full cohort. 2019 data cut. Medtronic data on file.
- ⁸ Melas N, Perdikides T, Saratzis A, Saratzis N, Kiskinis D, Deaton DH. Helical endostaples enhances endograft fixation in experimental model using human cadaveric aortas. *J Vasc Surg.* 2012; 55(6):1726-33
- ⁹ Chaudhuri A, Kim HK, Valdivia AR. Improved Midterm Outcomes Using Standard Devices and EndoAnchors for Endovascular Repair of Abdominal Aortic Aneurysms with Hyperangulated Necks [published online ahead of print, 2020 May 8]. *Cardiovasc Intervent Radiol.* 2020;10.1007/s00270-020-02488-4
- ¹⁰ Schlosser FJV, de Vries JPPM, Chaudhuri A. Is it time to insert endoanchors into routine EVAR? *Eur J Endovasc Surg.* 2017 Apr;53(4):458-459
- ¹¹ Endovascular Today Spotlight Article - EndoAnchors advance treatment of aortic aneurysms: Experience in a community hospital. Dr. Mazin Foteh. June 2016
- ¹² Litwinski RA, Donayre CE, Chow SL, Song TK, Kopchok G, Walot I, White RA., The role of aortic neck dilation and elongation in the etiology of stent graft migration after endovascular abdominal aortic aneurysm repair with a passive fixation device. *J Vasc Surg.* 2006 Dec;44(6):1176-81
- ¹³ Veith 2019 Dr. Jordan ANCHOR 4 year presentation: <https://www.veithsymposium.org/abstracts/2019/vei/739.pdf>
- ¹⁴ Stather PW, Sayers RD, Cheah A, Wild JB, Bown MJ, Choke E. Outcomes of endovascular aneurysm repair in patients with hostile neck anatomy. *EJVES.* 2012;44(6):556-561
- ¹⁵ Medtronic Presents 1-Year Valiant Navion and 4-Year ANCHOR Registry Outcomes. 22 November 2019. Endovascular Today. <https://evtoday.com/news/medtronic-presents-1-year-valiant-navion-and-4-year-anchor-registry-outcomes>

HELI-FX™ & HELI-FX™ THORACIC ENDOANCHOR™ SYSTEMS

Indications for Use

The Heli-FX™ EndoAnchor™ system is intended to provide fixation and sealing between endovascular aortic grafts and the native artery. The Heli-FX™ EndoAnchor™ system is indicated for use in patients whose endovascular grafts have exhibited migration or endoleak, or are at risk of such complications, in whom augmented radial fixation and/or sealing is required to regain or maintain adequate aneurysm exclusion. The EndoAnchor™ implant may be implanted at the time of the initial endograft placement, or during a secondary (i.e. repair) procedure.

Contraindications

Treatment with the Heli-FX™ EndoAnchor™ system is contraindicated for use in the following circumstances:

- In patients with known allergies to the EndoAnchor™ implant material (MP35N-LT)
- In conjunction with the Endologix Powerlink™* endograft

Warnings

- The long-term performance of the EndoAnchor™ implant has not been established. All patients should be advised endovascular aneurysm treatment requires long-term, regular follow-up visits to assess the patient's health status and endograft performance. The EndoAnchor™ implant does not reduce this requirement.
- The EndoAnchor™ implant and the Heli-FX™ EndoAnchor™ system have been evaluated via in vitro testing and determined to be compatible with the Cook Zenith™*, Cook Zenith™* TX2™*, Gore Excluder™*, Gore TAG™*, Medtronic AneuRx™, Medtronic Endurant™, Medtronic Talent™ AAA, Medtronic Talent™ TAA, Medtronic Valiant Xcelerant™, Medtronic Valiant™ Captivia™, and Medtronic Valiant Navion™ endografts. Use with endografts other than those listed above has not been evaluated.
- The performance of the EndoAnchor™ implant has not been evaluated for securing multiple endograft components together. Not securing EndoAnchor™ implants into aortic tissue could result in graft fabric damage, component separation, and resultant Type III endoleaks.
- The performance of the EndoAnchor™ implant has not been evaluated in vessels other than the aorta. Use of the EndoAnchor™ implant to secure endografts to other vessels may result in adverse patient consequences such as vascular perforation, bleeding, or damage to adjacent structures.
- The performance of the EndoAnchor™ implant has not been evaluated for securing multiple anatomical structures together. Such use could result in adverse patient consequences such as vascular perforation, bleeding, or embolic events.

MRI Safety and Compatibility

- The EndoAnchor™ implants have been determined to be MR Conditional at 3T or less when the scanner is in Normal Operating Mode with whole-body-averaged SAR of 2 W/kg, or in First Level Controlled Mode with a maximum whole-body-averaged SAR of 4 W/kg.
- Please refer to documentation provided by the endograft system manufacturer for MR safety status of the endograft system with which the EndoAnchor™ implants are being used.

Potential Adverse Events

Possible adverse events that are associated with the Heli-FX™ EndoAnchor™ system, include, but are not limited to:

- Aneurysm rupture
- Death
- EndoAnchor™ implant embolization
- Endoleaks (Type III)
- Enteric fistula
- Failure to correct/prevent Type I endoleak
- Failure to prevent endograft migration
- Infection
- Renal complications (renal artery occlusion/dissection or contrast-induced acute kidney injury)
- Stroke
- Surgical conversion to open repair
- Vascular access complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vessel damage, including dissection, perforation, and spasm

Please reference product *Instructions for Use* for more information regarding indications, warnings, precautions, contraindications and adverse events. Additional potential adverse events may be associated with endovascular aneurysm repair in general. Refer to the *Instructions for Use* provided with the endograft for additional potential adverse events.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a licensed healthcare practitioner. See package inserts for full product information.

CAUTION: EndoAnchor™ implant locations should be based upon a detailed examination of the preoperative CT imaging in cases involving irregular or eccentric plaque in the intended sealing zones. EndoAnchor™ implants should be implanted only into areas of aortic tissue free of calcified plaque or thrombus, or where such pathology is diffuse and less than 2mm in thickness. Attempting to place EndoAnchor™ implants into more severe plaque or thrombus may be associated with implantation difficulty and suboptimal endograft fixation and/or sealing.

[medtronic.com/aortic](https://www.medtronic.com/aortic)

UC202103224 EN ©2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. For distribution in the USA only. 08/2020

Medtronic