



Go Beyond Surveillance Biopsies with HeartCare

**The first multimodal test with studies demonstrating
a significant reduction in biopsies, resulting in more
patient-centric clinical decision making**




Sam D., heart transplant
recipient, and his wife

HeartCare Offers a **Comprehensive Approach** to **Post-transplant Surveillance**



A non-invasive surveillance tool with the dual power of cell-free DNA and gene expression profiling

- Measures mRNA of the recipient
- Measures immune response
- Originates from WBCs




- Measures dd-cfDNA from the transplanted heart
- Originates from the graft


HeartCare Has Proven Validation and Clinical Utility for Heart Transplant Rejection Surveillance



Now published: the inaugural publication from SHORE in JHLT²




Patients surveilled with HeartCare in SHORE had excellent outcomes²




- At 2-Years Post Transplant:**
- ✓ **94.9%** survival rate
 - ✓ **97.3%** of surviving patients had normal graft function
 - ✓ **61%** mean LVEF for surviving patients

Dual Positive HeartCare Results **Better Identify ACR** than A Single Test Alone²



Dual Positive HeartCare Results **Increased the Odds of ACR** by ~4X²





Positive Likelihood Ratio gives the change in the odds of having a diagnosis in patients with a positive test. A LR+ is mathematically defined as sensitivity / (1-specificity)

High AlloMap is: ≥ 30 for 2-6 months or ≥ 34 for >6 months | High AlloSure is $\geq 0.20\%$

HeartCare Interpretation Improved Over Time, Resulting in **More Patient-Centric Biopsy Decision Making**

Patients monitored with HeartCare had fewer biopsies over time, with fewest biopsies occurring in patients without dual positive HeartCare results.



Clinical Interpretation for HeartCare - ACR Surveillance²

HeartCare helps clinicians precisely identify patients who are at higher risk of ACR

HIGH ALLOMAP / LOW ALLOSURE A biopsy is unlikely to reveal ACR (ACR positivity in SHORE ² = 1.9%) Consider biopsy or repeat HeartCare testing earlier if: <ul style="list-style-type: none">AlloSure level is close to threshold and increased from prior measurement Consider other pathological causes of an increased AlloMap: <ul style="list-style-type: none">CMV infection <div>Medication adherence review recommended</div>	DUAL POSITIVE HEARTCARE (High AlloSure / High AlloMap) A biopsy is more likely to reveal ACR (ACR positivity in SHORE ² = 9.2%) <ul style="list-style-type: none">Biopsy should be considered <div>Medication adherence review recommended</div>
DUAL NEGATIVE HEARTCARE (Low AlloSure / Low AlloMap) A biopsy is unlikely to reveal ACR (ACR positivity in SHORE ² = 1.5%) <ul style="list-style-type: none">Continue protocol immuno-optimization	LOW ALLOMAP / HIGH ALLOSURE A biopsy is unlikely to reveal ACR (ACR positivity in SHORE ² = 4.3%) Consider a biopsy or repeat HeartCare testing earlier if: <ul style="list-style-type: none">AlloMap is close to threshold and AlloSure has increased by ≥0.20% from prior measurementRecent treatment for rejection (<21 days) or current prednisone >20 mgAt risk of Antibody Mediated Rejection/markedly elevated AlloSure Consider other possible pathological causes of an increased AlloSure: <ul style="list-style-type: none">Cardiac allograft vasculopathySevere infectionAntibody Mediated Rejection (AMR) / Donor specific antibodies <div>Medication adherence review recommended</div>


The table is provided for informational purposes only and is not intended as medical advice. A physician's test selection and interpretation, diagnosis, and patient management decisions should be based on his/her education, clinical expertise, current guidelines, and assessment of the patient. Please refer to publications for detailed clinical discussion. Clinical interpretation of AlloSure in Heart Transplantation graph is based on data from all commercial samples.

This table is designed for the context of surveillance testing for ACR; For patients that are at risk of AMR or being tested in other clinical context, different guidance may apply.

High AlloMap is ≥30 for 2-6 months or ≥34 for >6 months | High AlloSure is ≥0.20%


HeartCare Sample Report Offers an **Integrated View** to **Drive Informed Decisions**

- Individual test results (AlloMap or AlloSure Heart) are released when available on individual reports
- Combined HeartCare results are provided on an integrated report that includes both AlloMap and AlloSure Heart when both test results are available




HeartCare Helps Patient Management **Without the Limitations** of Traditional Biopsies

Heart transplant surveillance is complex



Traditional biopsies have significant limitations for rejection surveillance


- Resource intensive
- Risk of exposure to pathogens
- Sampling errors and limitations
- Invasive, painful, and frequent (sometimes >10 Biopsies in the 1st year post-transplant)
- Interobserver variability in interpretation




Almost 40% of rejections called by local pathologists **were downgraded** when reviewed by a panel of central pathologists

Agreement between local and central pathology for rejection ¹		
Local Grades	Central Grades	
	2R Rejection	< 2R Rejection
2R Rejection	46/76 (60.5%)	30/76 (39.5%)
< 2R	63/862 (7.3%)	799/862 (92.7%)

Traditional tools for surveillance monitoring may not allow for early intervention



2023 ISHLT guidelines support use of HeartCare in routine monitoring of heart transplant patients



References to ISHLT are offered solely to support AlloMap's FDA indications and should not be construed as supporting any other use. AlloMap should be solely used in conjunction with standard clinical assessment

A Suite of **Innovative Transplant Management Solutions**, Enabling More Cohesive Care

CareDx is dedicated to patients and providers throughout the transplant journey

The transplant journey is complex. Centers and providers are on a constant quest to deliver better outcomes while navigating the reality of operational challenges and fragmentation along the way.

CareDx provides a connected and innovative set of transplant management solutions to enable more cohesive care.



Patient Management

Solutions enabling providers to engage patients in their own care at every stage in the transplant journey



Specialty pharmacy dedicated to the unique needs of pre/post-transplant patients



Medication management, education and adherence platform



Free pre-to-post transplant app that provides medication management and biometric tracking to patients and is fully integrated with AlloHome and TxAccess



Customizable remote patient monitoring platform and monitoring service



Operational Support

Administrative and workflow support that enable process excellence and reduce inefficiencies and errors



Customizable administrative service



Referral management and communication tool



Quality Improvement

Advanced algorithms, predictive analytics and quality program tracking to help improve patient outcomes



Comprehensive transplant quality management software



Informed Decision Support

Data-driven insights that help providers assess individual patient's status and deliver tailored care



AiCAV helps to stratify patients based on their risk of developing CAV, a leading cause of heart transplant failure [In development]



Laboratory Management

Software and interoperability solutions for the histocompatibility and immunogenetics community




Laboratory Information Management System (LIMS) purpose-built for the world of histocompatibility and immunogenetics

References

1. Henriksen EJ, Moayed Y, Purewal S, et al. Combining donor derived cell free DNA and gene expression profiling for non-invasive surveillance after heart transplantation [published online ahead of print, 2022 May 12]. *Clin Transplant*. 2022;e14699. doi:10.1111/ctr.14699.
2. Khush KK, Hall S, Kao A, et. al. Surveillance with Dual Non-invasive Testing for Acute Cellular Rejection After Heart Transplantation: Outcomes from the Surveillance HeartCare Outcomes Registry (SHORE). *J Heart Lung Transplant*. 2024. <https://doi.org/10.1016/j.healun.2024.05.003>
3. Pham MX, Teuteberg JJ, Kfoury AG, et al. Gene-expression profiling for rejection surveillance after cardiac transplantation. *N Engl J Med*. 2010;362(20):1890-1900. doi:10.1056/NEJMoa0912965.
4. Kobashigawa J, Patel J, Azarbal B, et al. Randomized pilot trial of gene expression profiling versus heart biopsy in the first year after heart transplant: early invasive monitoring attenuation through gene expression trial. *Circ Heart Fail*. 2015;8(3):557-564. doi:10.1161/CIRCHEARTFAILURE.114.001658.
5. Khush KK, Patel J, Pinney S, et al. Noninvasive detection of graft injury after heart transplant using donor-derived cell-free DNA: A prospective multicenter study. *Am J Transplant*. 2019;19(10):2889-2899. doi:10.1111/ajt.15339. ; D-OAR is a sub-study of the Outcomes AlloMap Registry (OAR).
6. Velleca A, Shullo MA, Dhital K, et al. The International Society for Heart and Lung Transplantation (ISHLT) guidelines for the care of heart transplant recipients. *J Heart Lung Transplant*. 2023;42(5):e1-e141. doi:10.1016/j.healun.2022.10.015

HeartCare Provides **Actionable Information** During Surveillance



HeartCare Has Been **Clinically Validated in Multi-Center Prospective Studies** Including >4,500 Patients

2006	Cardiac Allograft Rejection Gene Expression Observation (CARGO) First multicenter validation of GEP test to identify heart transplant patients at low risk mod/severe rejection	<i>AJT</i>	N = 222
2010	Invasive Monitoring Attenuation through Gene Expression (IMAGE) Prospective, randomized, controlled, multi-center trial where Surveillance for rejection with AlloMap (gene expression profiling) was compared to routine biopsy (standard of care)	<i>NEJM</i>	N = 602
2015	Early Invasive Monitoring Attenuation through Gene Expression (eIMAGE) Randomized controlled trial comparing GEP to EMB as early as 55 days post transplant	<i>Circ: Heart Fail.</i>	N = 60
2016	CARGO II The CARGO II study validated the CARGO I study findings of GEP score performance in a predominantly European based cohort	<i>Eur Heart J.</i>	N = 472
2019	Outcomes AlloMap® Registry (OAR) Registry 5-year follow-up of patients from 2013-2019	<i>JHLT</i>	N = 1,504
2019	Donor-derived cell free DNA Outcomes AlloMap Registry (D-OAR) Subset of patients co-enrolled in OAR Prospective, multi-center registry with 26 centers to determine whether dd-cfDNA level in heart transplant recipients' blood can differentiate rejection from no rejection	<i>AJT</i>	N = 740
2024	Surveillance HeartCare Outcomes Registry (SHORE) Prospective registry captured events in patients during their first 5 years post-transplant	<i>JHLT</i>	N = 2732