

A close-up, black and white photograph of a metal packing coil, showing the intricate, tightly packed spiral of the wire. The coil is the central focus, with a blurred background that suggests it is part of a larger mechanical assembly. The lighting highlights the metallic texture and the repetitive pattern of the coil.

***swift*PAC**TM

Packing Coil

swiftPAC™

Designed for **Efficient, Dense Packing**



**.0135" Primary Diameter
Coil for Volume + Softness**

**Unique 2D
Sinusoidal Shape**

Designed to help the coil break
effortlessly and seek empty space

**Compatible with
10-System
Microcatheters**

10

60 cm

3 cm

Range of
Lengths from
**3 cm up to
60 cm**


**Unique coil shape
engineered for packing**

Photographs taken by and on file at Penumbra, Inc. One SwiftPAC 60 cm coil deployed with two 10-System anchoring coils. Bench test results may not be indicative of clinical performance.

Versatility Across Clinical Applications

Dr. Alan North

St. Dominic Jackson Memorial Hospital, MS




MMA Embolization

6 cm SwiftPAC coil (1)

Dr. Mohsen Nouri

Jamaica Hospital Medical Center, NY




CCF Embolization

60 cm SwiftPAC coil (2), 15 cm SwiftPAC coil (2)

Very Small Aneurysm

2.5 mm × 3 cm SMART COIL® Extra Soft + 3 cm SwiftPAC coil




Dr. Giuseppe Ammirati

Sharp Grossmont Hospital, CA

Flow Diverter and Coiling

30 cm SwiftPAC coil (1)



Dr. Nam Yoon

Intermountain Medical Center, UT

Case images used with permission. Consents on file at Penumbra, Inc. Individual results may vary depending on patient-specific attributes and other factors.

Ordering Information

Access Portfolio

Catalog Number	Description	Working Length (cm)	OD (in. / mm)	ID (in.)	Wire Compatibility (in.)
MIDWAY™ Intermediate Catheters					
43MWAY115	MIDWAY 43 Delivery Catheter, 115 cm	115	.060 / 1.52	.043	.038
43MWAY125	MIDWAY 43 Delivery Catheter, 125 cm	125	.060 / 1.52	.043	.038
62MWAY115	MIDWAY 62 Delivery Catheter, 115 cm	115	.076 / 1.93	.062	.038
BENCHMARK™ 071 Access System					
			(F)		
BMK6F95	BENCHMARK 071 Access Catheter, 95 cm	95	6 / 2.11	.071	.035 – .038
BMK6F95M	BENCHMARK 071 Access Catheter, 95 cm MP	95	6 / 2.11	.071	.035 – .038
BMK6F105	BENCHMARK 071 Access Catheter, 105 cm	105	6 / 2.11	.071	.035 – .038
BMK6F105M	BENCHMARK 071 Access Catheter, 105 cm MP	105	6 / 2.11	.071	.035 – .038
BENCHMARK BMX®81 Access System					
BMX7F95	BENCHMARK BMX81 Delivery Catheter, 95 cm	95	7 / 2.46	.081	.038
BMX7F105	BENCHMARK BMX81 Delivery Catheter, 105 cm	105	7 / 2.46	.081	.038
BMX7F115	BENCHMARK BMX81 Delivery Catheter, 115 cm	115	7 / 2.46	.081	.038

SwiftPAC Coil

Packing Coil



Catalog Number	Length (cm)
400SWPC03	3
400SWPC06	6
400SWPC10	10
400SWPC15	15
400SWPC30	30
400SWPC45	45
400SWPC60	60


Detachment Handle



Catalog Number	Description
SDH1	Swift™ Coil Detachment Handle

PAC400™ Coil


Pack Behind POD400 Backstop



Catalog Number	Description	Length (cm)
400PAC7	PAC400 Coil 7 cm	7
400PAC10	PAC400 Coil 10 cm	10
400PAC15	PAC400 Coil 15 cm	15
400PAC25	PAC400 Coil 25 cm	25
400PAC35	PAC400 Coil 35 cm	35
400PAC40	PAC400 Coil 40 cm	40
400PAC50	PAC400 Coil 50 cm	50
400PAC60	PAC400 Coil 60 cm	60

POD400™ Coil

Vessel Sacrifice



Catalog Number	Description	Target Vessel (mm)	Length (cm)
400POD3	POD400 Coil 3 mm	3	20
400POD4	POD400 Coil 4 mm	3 – 4	30
400POD5	POD400 Coil 5 mm	4 – 5	30
400POD6	POD400 Coil 6 mm	5 – 6	50



For the complete Penumbra IFU Summary Statements, please scan QR code or visit: peninc.info/risk

Photographs taken by and on file at Penumbra, Inc. Case images used with permission. Consents on file at Penumbra, Inc. Individual results may vary depending on patient-specific attributes and other factors.

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use (IFU) for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2024 Penumbra, Inc. All rights reserved. The Penumbra P logo, Swift, SwiftPAC, Penumbra SMART COIL, MIDWAY, BENCHMARK, BMX, POD, POD400, and PAC400 are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. 29091, Rev. A 08/24 USA

Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

penumbrainc.com