




Every Signal. Every Second.



VitalPatch®



Built for Professionals. Made for Patients.¹



Thorough Detection. Effortless Wear.



Our unique FDA-cleared technology helps accelerate informed clinical decision-making from rhythm to resolution^{2,3,4} while also delivering a simplistic and comfortable experience for patients.

Powered by our proprietary FDA-cleared algorithm that detects 4 key vital signs⁴:

- Heart Rate
- Heart Rate Variability
- Respiratory Rate
- Activity

One Patch. Multiple Modalities.



Our single-use biosensor helps clinicians deliver better cardiac care with the detection of key vital signs comfortably and continuously, in any setting² for up to 30 days.^{2,3}



4-in-1 cardiac monitor

Continuously captures patient data for multiple modalities in one device — including Mobile Cardiac Telemetry (MCT), Extended Holter and Event.



User-friendly experience and comfortable for patients

The lightweight, disposable patch eliminates hassle — so patients can wear and go about life with confidence.



FDA-cleared & performance tested^{1,2}

Wearable technology that assists physicians confidently care for their patients.



No Abrading Needed¹



Shower & Exercise Normally⁵

55.OB+

Total Number
of Heart Beats¹

220.OK+

Arrhythmia
Detected⁶



The Architecture of VitalPatch¹

ECG Performance	Continuous ECG, with cloud-based algorithm detection of 21 unique arrhythmias
Vital Signs	Heart Rate (HR) Heart Rate Variability (HRV) Respiratory Rate (RR) Activity (Steps) Activity (Movement, mG)
Weight	13g (0.46 ounces)
Dimensions	120mm x 41mm x 9.5mm (4.7" x 1.6" x 0.38")
Adhesive	Hypoallergenic hydrocolloid, latex-free.
Battery	Zinc-air, 168-hour life (7 days), no charging required.
Water Resistance	Patients have the freedom to go about their normal daily activities, such as showering. Swimming and submersion are not supported.
Operating	50 to 104°F (10 to 40°C), 10 to 95% relative humidity
Storage	32 to 104°F (0 to 40°C), 10 to 95% relative humidity
Shelf Life	9 months
Connectivity	Bluetooth Low-Energy (BLE), 10 feet line-of-sight between biosensor and phone.

CPT Codes⁷

Modality	Global	Hook-Up	Technical	Professional
Mobile Cardiac Telemetry (MCT)			93229	93228
3-7 days Holter	93241	93242	93243	93244
7+ days Holter	93245	93246	93247	93248
Cardiac Event Monitoring (CEM)	93268	93270	93271	93272

1.Data on file. VitalConnect 2025.

2.VitalConnect solutions are intended to be used for outpatient cardiac telemetry and in-patient monitoring in non-critical healthcare settings for non-urgent clinical decision-making.

3.May be affected by data transmission. Potential transfer data interruptions are listed in the IFU.

4.In the detection of key vital signs with Live Look-In™ ECG.

5.VitalPatch is water-resistant not waterproof. No swimming with the patch or submerging it in water.

6.Per Physician's MCT/CEM Notification Criteria

7.Information contained in this publication is not to be construed as legal or billing advice. CPT® is a registered trademark of the American Medical Association. All CPT® information provided in this publication is intended for illustration purposes only, and should be independently verified prior to billing application.