

ECGs WITH CLINICAL INSIGHTS

VitalPatch™

Mobile Cardiac Telemetry (MCT) Service

VitalConnect®

AN iRHYTHM COMPANY

Built for Professionals


Stay Connected to What's Vital

Secure Live Look-In*

Access a patient's full disclosure and data during or after their wear period as alerts arise or when symptoms are reported.


Historical Trends

As alerts arise or when symptoms are reported, providers can review historical trends.



VistaCenter® with Live Look-In*

FULL DISCLOSURE





LIVE LOOK-IN

*Live Look-In data may take 1–2 minutes to display after logging in to VistaCenter™ and may take up to 5 minutes. It does not provide real-time or near-real time ECG data. Refer to the VistaCenter Instructions for Use for additional information.


Get the Answers You Ordered

Representative ECG strips: For the fastest arrhythmia (AFib, VT, SVT) and slowest or clinically actionable (Pause, Block, or AIVR/IVR).

Onset, Duration, and Offset: Strips provided for all arrhythmia events.

Patient-Reported Symptoms: Detailed view of patient-reported symptoms with representative strips.

Ask your local representative for a sample end-of-service report.†



†End-of-service report processed by a third-party algorithm.

NOT ALL MONITORING
SERVICES ARE THE SAME

Choosing the Right Modality

VitalConnect provides four monitoring modalities through a single, lightweight 13 gram VitalPatch® monitor*, enabling clinical flexibility without the need to maintain multiple devices.

The VitalPatch® monitor is indicated for continuous collection of physiological data including the ECG data for up to 7 days.



*VitalPatch is replaced every 7 days during the prescribed monitoring period. Refer to the VitalPatch Instructions for Use for additional information.

Mobile Cardiac Telemetry (MCT)

SERVICE PRESCRIBED UP TO 30 DAYS*

- Live telemetry analyzed continuously.
- Events verified by qualified cardiac technicians.
- Actionable notifications relayed during wear.
- Continuous ECG analysis on end-of service report.†

Cardiac Event Monitor (CEM)‡

SERVICE PRESCRIBED UP TO 30 DAYS*

- Episodic reporting based on auto-triggered detection or patient-triggered symptom events.
- Event verification by qualified cardiac technicians and relayed during wear.
- Does not offer full ECG analysis on end-of-service report.†

Extended Holter (EH)

SERVICE PRESCRIBED UP TO 15 DAYS*

- Continuous ECG recording with retrospective analysis.
- No live telemetry or notifications during wear.
- Patch change after 7 days.*

24–48 Hour Holter

DURATION: 1–2 DAYS§

- Retrospective ECG monitoring for up to 48 hours.
- No live telemetry or notifications during wear.

INSURANCE FLEXIBILITY

Coverage That Works for You

The VitalPatch® monitor supports MCT and CEM services, providing flexibility based on patient coverage by enabling a physician-ordered transition from MCT to CEM without requiring the patient to change monitor.*†‡

*VitalPatch is replaced every 7 days during the prescribed monitoring period. Refer to the VitalPatch Instructions for Use for additional information. | ‡Services may be prescribed for up to 30 days.
†The VitalPatch monitor is indicated for continuous collection of physiological data including the ECG data for up to 7 days. | §Available upon request.

Made for Patients



PATIENT COMFORT

VitalPatch® monitor uses a hydrocolloid adhesive designed for patient comfort.

SERVICE UP TO 30 DAYS

Supports MCT service prescribed for up to 30 days. No patch charging required. Just replace the VitalPatch® monitor every 7 days.*

MULTIPLE LANGUAGES

Patient support resources in English and Spanish through the VistaPoint™ app on the VistaPhone™ relay device.

NOW GREEN!

The new, reusable electronics module reduces e-waste when returned.

DESIGNED FOR PATIENT COMFORT & FREEDOM

VitalPatch® monitor is discreet and seamlessly integrates into daily routines.[1]

13 GRAMS

Lightweight for patient comfort.[1]

DAILY LIFE

Patients can continue with their regular routines including showering and exercising.†

*The VitalPatch monitor is indicated for continuous collection of physiological data including the ECG data for up to 7 days. | †Minimize exposure directly under the shower head, excessive contact with soap, or scrubbing. Gently dry the device after showering. Do not submerge the device or use in a sauna. Please refer to the VitalPatch Instructions for Use for the full set of details. | [1]VistaPoint Instructions for Use. VitalConnect.

Technology at the Core

Behind the Monitoring Service

Numbers Behind the Algorithm

The VitalRhythm™ Algorithm[2] was trained on a database of over 23,000 patients from 354 US clinical sites.[3]



Advanced Ambulatory Cardiac Monitoring Service


The VitalPatch® monitor is supported by our ASIC chip's ability to synthesize 8+ physiologic parameters simultaneously.[1]

Seamless IT Integration

COMPATIBLE WITH MAJOR EHR SYSTEMS

Including Epic®, Oracle® (formerly Cerner®), Meditech®, and other top systems. This bidirectional integration eliminates the need for manual data transfers or extra logins.

[1]VistaPoint Instructions for Use. VitalConnect. | [2]The VitalRhythm Algorithm is FDA cleared for during wear ECG data analysis | [3]Data on file. VitalConnect.



VitalConnect, Inc.
2870 Zanker Road, Suite 100 | San Jose, CA 95134
(408) 963-4600 | vitalconnect.com

Contact Us Today

VitalConnect®
AN iRHYTHM COMPANY

© 2026 VitalConnect. A list of our trademarks can be found at www.vitalconnect.com.
All other marks may be trademarks of their respective companies. | MAT-0458-1